Instructions for Use

Instructions for Use – GB
Gebrauchsanweisung – DE
Instructions d'utilisation – FR
Bruksanvisning – SE
Istruzioni per l'uso – IT
Instrucções para a utilização – PT
Οδηγίες Χρήσης – GR
Manual de instrucciones – ES
Brugsanvisning – DK
Käyttöohjeet – FI
Novadila za uporabo – SL

ABRIS ABR screening module *the Eclipse black box*

Valid from software version 1.04
80799504 – 10/2010

Instructions for Use GB

ABRIS ABR screening module

- for the Eclipse black box

Intended Use

The ABRIS ABR infant screening software module allows for a quicker ABR testing time than traditional ABR and therefore a more efficient Universal Hearing Screening program. However, an immediate and comprehensive OAE or ABR follow up can be made within minutes on babies failing the initial screen, if you also have the EP15/EP25 ABR system and/or the TEOAE25 installed on your Interacoustics Eclipse black box.

The measurement can be performed using traditional electrodes and insert phones or TDH39.

Precautions

| | |
|--|--|
|  | <p>WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.</p> |
|  | <p>CAUTION, used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.</p> |
|  | <p>NOTICE is used to address practices not related to personal injury.</p> |

1. Be sure to use only stimulation intensities, which will be acceptable for the patient.
2. The transducers (headphones) supplied with the instrument are calibrated to this instrument - exchange of transducers requires a new calibration.
3. It is recommended that parts which are in direct contact with the patient (e.g. earphone cushions) are subjected to standard disinfecting procedure between patients. This includes physically cleaning and use of a recognized disinfectant. Individual manufacturer's instruction should be followed for use of this disinfecting agent to provide an appropriated level of cleanliness.
4. Please note that the CE Marking is only legal if the CE Instruction is translated into the national language of the user no later than at the delivery to him, if the national legislation demands a text in the national language according to MDD article 4.4.
5. Although the instrument fulfils the relevant EMC requirements precautions should be taken to avoid unnecessary exposure to electromagnetic fields, e.g. from mobile phones etc. If the device is used adjacent to other equipment it must be observed that no mutual disturbance appears.
6. Please note that if connection is made to standard equipment like printers and network, special precautions must be taken in order to maintain medical safety. Optical isolation unit for USB is available from your supplier.

7. The specification for the instrument is valid if the instrument is operated within the following environmental limits:

Temperature: 15°C to 35 °C.

Humidity: 30 %RH to 90 %RH


Supply voltage: 100 – 240 Vac

Supply frequency: 50 – 60 Hz

If this apparatus is connected to one or more other devices with medical CE marking, to make up a system or pack, the CE marking is only valid also for the combination if the supplier has issued a declaration stating that the requirements in the Medical Device Directive article 12 are fulfilled for the combination.

Making a recording using electrodes

- 1) To enter test mode from the OtoAccess™, double click the ABRIS icon 
- 2) Apply electrodes. Blue and Red for Left and Right mastoid respectively. White for vertex (or hairline at forehead) and Black for lower forehead. Most neonates with clean (but not oiled) skin do not need abrasive cleaning prior to electrode montage.
- 3) Check Impedance with PreAmplifier in "Imp" mode. When turning the knob counter clockwise, the lights should go green somewhere below the 3k Ω indication.
- 4) Otherwise you may want to redo electrode montage (and maybe clean or abrase skin) at electrode site in question.
- 5) Click on Imp botton again to return to ERA-mode.
- 6) Observe EEG level – too much EEG noise from an unrest full baby (or from a very poor electrode montage) will cause noise bar to be in the red zone (or the EEG to be red, if the raw EEG is displayed instead of the noise bar).
- 7) Once noise levels are OK, press START.
- 8) Wait for test to complete.
- 9) Save and Exit by pressing .

Connection Panel Eclipse

| Position: | Symbol: | Function: |
|-----------|----------------|---|
| 1 | Power | Turns power on/off |
| 2 | Mains 50-60 Hz | Socket for mains cable |
| 3 | AUX | For future use |
| 4 | Preamp. | Socket for preamplifier |
| 5 | Pat. Resp. | Socket for patient response knob |
| 6 | Trigger In/Out | Socket for trigger input/output |
| 7 | Talk Back | Socket for talk back microphone |
| 8 | Talk Forward | Socket for talk forward microphone |
| 9 | OAE | Socket for the OAE probe |
| 10 | Left | Socket for left phone |
| 11 | Right | Socket for right phone |
| 12 | Bone | Socket for bone conductor |
| 13 | USB/PC | Socket for USB cable or PC |
| 14 | DC | Socket for power supply for optical USB extension cable |
| 15 | | Connection for functional ground |

Explanation to symbols which can be found on the instrument:

| | |
|---|---|
| I | On (Power: connection to the mains) |
| O | Off (Power: disconnection from the mains) |
| ⊥ | Ground |
| ⚠ | See explanation in manual |
| ⚡ | Type BF equipment |

Technical Specifications

ABRIS Specifications:

| | |
|-----------------------------------|--|
| Medical CE-mark: | The CE-mark indicates that Interacoustics A/S meets the requirements of Annex II of the Medical Device Directive 93/42/EEC. Approval of the quality system is made by TÜV – identification no. 0123. |
| System Platform: | The same as EP15/25 ABR system and TEOAE25, DPOAE20 screening and diagnostic OAE |
| Hardware Platform: | ECLIPSE black box with USB connection to laptop |
| Database: | OtoAccess™ |
| Stimulus: | Click 93Hz stimulus rate |
| Level: | 30dBHL, 35dBHL or 40dBHL |
| Security: | Password protection of test parameters |
| Test time: | Typically less than 20 seconds per ear |
| Algorithmical Sensitivity: | 99.99% |
| Specificity: | Typically 97% |

Eclipse Specifications only:

| | |
|---------------------------------|--|
| Power Supply: | Input Volts: 100 – 120V, 200 – 240V, 50/60 Hz |
| | Consumption: 60VA |
| | Fuses: T 2A L |
| | Safety: IEC 60601-1 |
| Minimum PC Requirements: | 256MB RAM, 1 GB harddisk, Windows98se, WindowsXP or Windows2000, Vista USB 1.1 or better, Pentium III 800 MHz. |
| Construction: | Metal cabinet |
| Dimensions: | (L x W x H) 28 x 32 x 5,5 cm |
| Weight: | 2,5 kg / 5,5 lbs excluding accessories |

Detachable Parts

ABRIS:

- ABRIS Software
- Operation Manual, CE Manual
- Pre-amplifier and EARTone ABR Insert Earphones
- Eclipse as well as OtoAccess™ software (if the software is not ordered as a kit for upgrading an existing EP15/25 or TEOAE25, DPOAE20 system).

Concerning Repair and technical assistance

1. Interacoustics is only considered to be responsible for the validity of the CE marking, effects on safety, reliability and performance of the equipment if:

assembly operations, extensions, readjustments, modifications or repairs are carried out by authorised persons,

a 1 year service interval is maintained

the electrical installation of the relevant room complies with the appropriate requirements, and

the equipment is used by authorised personnel in accordance with the documentation supplied by Interacoustics.

2. It is important that the customer (agent) fills out the RETURN REPORT, located in the back of the operational manual every time a problem arises and sends it to Interacoustics A/S, Drejervaenget 8, DK-5610 Assens, Denmark. This should also be done every time an instrument is returned to Interacoustics. (This of course also applies in the unthinkable worst case of death or serious deterioration to patient or user).

3. When instrument fuses need renewal, the correct type as stated on the instrument shall be used.

General maintenance procedure

The performance and safety of the instrument will be kept if the following recommendations for care and maintenance are observed:

It is recommended to let the instrument go through at least one annual evaluation, to ensure that the acoustical, electrical and mechanical properties are correct. This should be made by an experienced workshop in order to guarantee proper service and repair.

Before the connection to the mains, be sure that the local mains voltage corresponds to the voltage labelled on the instrument.

Observe that no damage is present on the insulation of the mains cable or the connectors and that it is not exposed to any kind of mechanical load, which could involve damage.

For maximum electrical safety, turn off the power from a mains powered instrument when it is left unused.

Do not place the instrument next to a heat source of any kind, and allow sufficient space around the instrument to ensure proper ventilation.

To ensure the reliability of the instrument, periodic biological measurements should be performed on a person with known data. This person could be the operator him/herself.

If the surface of the instrument or parts of it is dirty, it can be cleaned using a soft cloth moistened with a mild solution of water and dish washing cleaner or similar. The use of organic solvents and aromatic oils must be avoided. Always disconnect the mains plug during the cleaning process, and be careful that no fluid is entering the inside of the instrument or the accessories.

After each examination of a patient, proper cleaning must ensure that there is no contamination on the parts in connection with patients. General precautions must be observed to prevent transmission of disease from one patient to another. If the ear cushions or ear tips are contaminated, it is strongly recommended to remove them from the transducer before they are cleaned. Frequent cleaning using water may be used, but periodic use of a mild disinfectant may also be used. The use of organic solvents and aromatic oils must be avoided.

Great care should be exercised by the handling of earphones and other transducers, as mechanical shock may cause change of calibration.

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