Instructions for Use

Instructions for Use - GB
Gebrauchsanweisung - DE
Instructions d'utilisation - FR
Istruzioni per l'uso - IT
Manual de instrucciones - ES
Instruções de utilização PT
Gebruiksaanwijzing - NL
Brugsanvisning - DK

Rotary Chair Nydiag 200

80705602 ver. 05/2009

Instructions for Use - GB

Rotary Chair Nydiag 200

General information

The Rotary Chair is constructed to a maximum patient weight of 135 kilos

This instruction manual is a part of the instrument and should always be kept at hand.

It documents the state of development of the unit, valid at the time of printing, and contains all the information required by the IEC 60601-1.

Exact adherence to the instruction manual is a prerequisite for the perfect and intended functioning of the Interacoustics instrument.

Always be sure to disconnect the power supply before cleaning or disinfecting mains operated instruments.

Interacoustics A/S only guarantees for the safety, reliability and function of the instrument, if: installation, extensions, modifications and repair are exclusively carried out by personnel authorized by Interacoustics A/S;

the room in which the equipment is installed complies with the country-specific regulation;

the installation site complies with the ambient conditions of the instrument;

the unit is used in accordance with the instruction manual after initiation in the operation;

The device is not approved for operation in potentially explosive rooms or flammable atmospheres.

If fluid gets into the instrument, the instrument will have to be switched off at once. The fluid must immediately be removed by Interacoustics service technicians and the instrument tested for technical safety. Defective connectors and leads must immediately be replaced or repaired by a specialist or by personnel authorized by Interacoustics A/S.

Precautions

Following the ANSI recommendations (American National Standards Institute) for safety notes, specific passages of this instruction manual are clearly marked as safety notes.

	WARNING indicates a hazardous situation which, if not avoided, could result in death or serious injury.
	CAUTION , used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury.
	NOTICE is used to address practices not related to personal injury

Hardware installation

The rotary chair must be secured to the floor with 4 screws through the bottom rim of the base. Please use 14 mm dowels and 10 mm screws.

When orienting the chair in the examination room please observe that the default zero position of the chair is indicated by the exit point of the cables from the chair.

Never connect the instrument to the power outlet before the side panels and the top cover have been correctly refitted and secured.

Consult with the product specification plate on the unit for compliance between the unit's own data and those of the local power supply system (mains voltage and mains frequency) before connecting the unit.

Connect only if all data agrees!

Inspect the mains connection cables for visible damage prior to establishing the connection. Damaged cables or plugs must immediately be replaced by Interacoustics service technicians, an authorized person or specialist.

Connection Panel

The backrest features a connection panel with purpose specific connections:

2 x 6-pM FireWire ports
1 x USB 2.0 port (not active)
1 x DC out (not active)

Software Installation

The VNG installation CD received with the VNG system holds the relevant software for the CanBus-to-USB adapter and the Rotary Chair Control software.

The Rotary Chair can be connected to any laptop or desktop PC with the following specifications:

- Intel Core 2 Duo processor 1.8 GHz or better
- Minimum 1 GB of RAM
- Laptop PC: One 34 mm PCExpressCard slot available
- Desktop PC: Texas Instruments chipset PCI FireWire® board with dual 6-pM ports
- Three USB ports available
- One VGA port available
- Monitor resolution 1024 x 768 or better
- WindowsXP™ Operating system

Starting the program

After switching on the power of the rotary chair, double-click the icon

on your Desktop. The following window appears.

Functions:

<input type="radio"/> Set Position	On the left hand side of the setup window you will find the five different modes of chair movements:
<input type="radio"/> Set Velocity	
<input type="radio"/> Sinusoidal	
<input checked="" type="radio"/> Triangular	
<input type="radio"/> Trapeze	

- Set position: Sets the chair to the angle defined in the setup
- Set velocity: Starts to move the chair to a certain velocity after pressing the start-button
- Sinusoidal: Sinusoidal smooth pursuit movement
- Triangular: Triangular smooth pursuit movement
- Trapeze: The chair velocity follows a trapeze pattern. This is used for velocity-step test and similar rotary tests.

1. Set Position Tab
2. Set Velocity Tab
3. Sinusoidal Tab
4. Triangular Tab
5. Trapeze Tab

Set position function

- In the **Set Position** mode the chair is moved to a pre-defined position; a default "Zero Position". After a test the chair will move to this position.
- If the examiner wants to define a second "Zero Position" where the chair faces a screen or similar, this position is defined in relation to the default Zero Position.
- Overview:

The **Position p** can be set between -180° ... 180° . **Velocity v** sets the velocity of the movement. **Acceleration a** sets the acceleration until **Velocity v** is reached and the deceleration from **Velocity v** until the chair stops.

- Comments:
 - a. After a test the chair returns to the default zero position which we define as zero degrees. As an example we now want the patient to face an object which is 180° to the right (a half cycle) of the zero position. The value '180' is entered in the window, and the 'Start' button is activated. The chair must perform a rightward half cycle rotation.

- b. For leftward movements towards the second zero position simply add a 'minus' sign before the value, like in "-180°"
- c. Velocity: low values such as 20°/s produce slow movements, whereas high values such as 100 °/s produce fast movements.

- **Testing the Set Position Function:**

- **Position – target of chair movement**
- **Velocity – speed of rotation to the intended position**
- **Acceleration – time to reach velocity**

First, press "Zero Position" button in the lower left corner of the Rotary Chair Control window.

1. **Position:** Enter a value such as 45, 90, 135 or similar in the **Position p:** box and press "Start".
Measure the number of degrees of movement – must equal the value entered in **'Position p:'** box.
2. **Velocity:** Check that velocity changes with different values entered in the **Velocity v:** box – low values produce slow movements and high values produce fast movements.
Enter the value '20' in the **'Velocity v:'** box and press "Start" button.
Allow the chair to accelerate up to set velocity (approx. 2 seconds with acceleration set at 10°/s²) and measure the time consumed for a full circle. Should be 18 seconds at velocity of 20 °/s
In order to obtain sufficiently long cycles it may be necessary to enter high values such as 720 in **Position p:** box – this will give at least 2 full circle rotations for measurement.
3. **Acceleration:** check that low values (such as 10°/s²) produce soft acceleration, whereas high values (such as 50°/s²) produce abrupt acceleration.
Enter the following values:

i. <i>Position</i>	=	720
ii. <i>Velocity</i>	=	100
iii. <i>Acceleration</i>	=	20

 Expected acceleration time is 5 seconds – check using a stop watch

NOTICE

Maximum acceleration is 100°/S² and must never exceed this value

Set velocity function

- In the **Set Velocity** mode the chair is rotated at a constant, pre-defined velocity until it is stopped via the "Stop" button.
- Once the acceleration phase is completed the chair will have reached its target velocity which it will keep until the "Stop" button is activated.
- Relevant VO25 test protocol: Data from **Set Velocity** can be recorded and analyzed in the StepRotation test
- VNG toolbar icon: 

- Overview:

The **Velocity v** can be set between -200°/s ... 200°/s. If the velocity is negative, the chair moves anti-clockwise, while at positive velocity-values it moves clockwise. **Acceleration a** sets the acceleration until this velocity is reached and the deceleration from the **Velocity v** until the chair stops.

- Comments:
 - A velocity of 90°/s means that the chair will do 1/4 or 0,25 full rotations per second – or take 4 seconds to do a full 360° cycle.
 - An acceleration value of 10°/s² indicates an acceleration time of 9 seconds before the set velocity of 90°/s is reached
- **Testing the Set Velocity Function:**
 - **Velocity °/s – constant velocity**
 - **Acceleration °/s² – linear increase in velocity up to constant velocity**
 1. A stopwatch is used to measure the velocity. Enter the value '90' in the **Velocity v:** box and press "Start". Allow the chair to accelerate to the set velocity and with the stopwatch, measure how long it takes for the chair to make a 360° cycle. At a constant speed of 90°/s it will take 4 seconds to do a full 360° cycle.
 2. Check that acceleration is fast with high values and moderate with low values

NOTICE

Maximum acceleration is 100°/S² and must never exceed this value

Set Sinusoidal test function (smooth pursuit, SHAT)

- In the **Sinusoidal** mode the system performs the movement which is used in the Smooth Pursuit test and in the SHAT test. The Smooth Pursuit test requires sinusoidal or triangular movements at a pre-defined frequency and amplitude.
- Purpose of **Sinusoidal**: The chair will rotate in pendular movements to the left and to the right or vice versa. Either with constant amplitude or with diminishing amplitude.
- Stimulus: Acceleration up to peak velocity, followed by deceleration of similar proportions.
- Relevant VNG test protocol: Data from **Sinusoidal** can be recorded and analyzed in the SinusPendular test
- VNG Toolbar icon: 

Conversion of Frequency to Acceleration

In order to avoid exceeding the maximum acceleration performance during higher frequencies of Slow Harmonics Acceleration Testing, please refer to the following conversion table for observation of upper limits:

Frequency	Acceleration @ Peak Velocity 60°/s	Acceleration @ Peak Velocity 50°/s	Acceleration @ Peak Velocity 40°/s
Hz	°/s ²	°/s ²	°/s ²
0.01	3.8	3.1	2.5
0.02	7.5	6.3	5.0
0.04	15.1	12.6	10.1
0.08	30.2	25.1	20.1
0.16	60.3	50.3	40.2
0.32	Not possible	100.5	80.4

SHAT test

For the SHAT test the frequency increases after each period. A commonly used procedure employs oscillation frequencies of 0.01, 0.02, 0.04, 0.08, 0.16 and 0.32 with peak angular velocities of 50°/sec at each frequency. The patient undergoes multiple cycles of oscillation at each frequency.

- **Testing the Sinusoidal Function:**

- **Velocity °/s – peak angular velocity between turning points**
- **Frequency Hz – duration of cycles between turning points/cycles**
- **Frequency change – steady or changing period length for each cycle**
- **Velocity change – steady or changing velocity for each cycle**
- **Number of cycles**

The various parameters of the sinusoidal test are difficult to test objectively. Apart from the number of cycles that can be counted, only subjective assessment is possible:

1. **Velocity:** does the chair move faster with increasing values for velocity? Compare values 20°/s with 50°/s – with Frequency Hz, Frequency change df and Velocity Change $dv = 1$
2. **Frequency:** With velocity fixed at a given value, does the period length increase with decreasing values for Frequency Hz? Compare with values '0,10' and '0,02' in the **Frequency Hz** box and press "Start" – at '0.10' the period should be shorter than at '0.02'
3. **Frequency Change:** values > 1 produce smaller cycles in terms of degrees rotated, whereas values < 1 produce longer cycles in terms of degrees rotated
4. **Velocity Change dv :** values > 1 produce faster movements whereas values < 1 produce slower movements
5. **Cycles:** check that the chair rotates the number of cycles ordered. One cycle is half the pendular movement; i.e. from zero, acceleration phase and then de-acceleration into zero again.

NOTICE

Maximum acceleration is 100°/S² and must never exceed this value

Hygiene

Always be sure to pull the plug before cleaning or disinfecting mains operated instruments.

The instrument must be cleaned with a non-fuzzy, slightly moistened piece of cloth (please do not soak it!) Only use a mild alcoholic fluid as a disinfectant.

When cleaning the instrument, any moisture, like e.g. condensation water, must necessarily be prevented from getting inside!

Chemicals required for operation or care of the unit must always be stored, prepared and made available in specially labelled containers to prevent any mistakes.

Maintenance

Before switching on the instrument make sure that the mains cable, mains plug, socket and input at the instrument are in a perfect state. If the drive of the rotary chair becomes noisy the driving belts should be renewed. We recommend renewing the driving belts after five years.

The operator is responsible for the perfect state of the instrument. It is recommended by Interacoustics to have safety checks carried out at one-year intervals.

Immediate maintenance is required if:

**the instrument was subjected to extreme mechanical stress (impact, defective cable due to inadmissible traction);
fluid gets into the instrument;
cables and/or connectors are defective;
rubber-joints have developed cracks.**

For further details please contact your Interacoustics representative.

SYMBOLS ON INTERACOUSTICS COMPONENTS

Ground conductor connection

Unit of type B

Technical Specifications

Maximum speed	200 deg/s
Maximum acceleration	100 deg/s ²
Maximum patient weight	135 kilos
Reclining backrest	manually operated from 0 deg (horizontal) to 90 deg (upright)
Emergency stop	emergency stop button disconnects the motor power
Patient alarm button	sends an alarm signal to the computer and turns off power to the motor
Weight	175 kg
Shipping weight	210 kg
Dimensions	90 cm x 70 cm x 160 cm
Shipping dimensions	100 cm x 120 cm x 190 cm
Power supply	110-230 V~ (50/60 Hz)/ 4A max

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Funzioni:

<input type="radio"/> Set Position <input type="radio"/> Set Velocity <input type="radio"/> Sinusoidal <input checked="" type="radio"/> Triangular <input type="radio"/> Trapeze	<p>On the left hand side of the setup window you will find the five different modes of chair movements:</p> <ul style="list-style-type: none"> - Set position: Sets the chair to the angle defined in the setup - Set velocity: Starts to move the chair to a certain velocity after pressing the start-button - Sinusoidal: Sinusoidal smooth pursuit movement - Triangular: Triangular smooth pursuit movement - Trapeze: The chair velocity follows a trapeze pattern. This is used for velocity-step test and similar rotary tests.
--	--

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Funciones:

<ul style="list-style-type: none"><input type="radio"/> Set Position<input type="radio"/> Set Velocity<input type="radio"/> Sinusoidal<input checked="" type="radio"/> Triangular<input type="radio"/> Trapeze	<p>On the left hand side of the setup window you will find the five different modes of chair movements:</p> <ul style="list-style-type: none">- Set position: Sets the chair to the angle defined in the setup- Set velocity: Starts to move the chair to a certain velocity after pressing the start-button- Sinusoidal: Sinusoidal smooth pursuit movement- Triangular: Triangular smooth pursuit movement- Trapeze: The chair velocity follows a trapeze pattern. This is used for velocity-step test and similar rotary tests.
--	--

1. Ficha Set Position (Establecer posición)
2. Ficha Set Velocity (Establecer velocidad)
3. Ficha Sinusoidal (Sinusoidal)
4. Ficha Triangular (Triangular)
5. Ficha Trapeze (Trapezio)

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Functies:

<input type="radio"/> Set Position	On the left hand side of the setup window you will find the five different modes of chair movements:
<input type="radio"/> Set Velocity	- Set position: Sets the chair to the angle defined in the setup
<input type="radio"/> Sinusoidal	- Set velocity: Starts to move the chair to a certain velocity after pressing the start-button
<input checked="" type="radio"/> Triangular	- Sinusoidal: Sinusoidal smooth pursuit movement
<input type="radio"/> Trapeze	- Triangular: Triangular smooth pursuit movement
	- Trapeze: The chair velocity follows a trapeze pattern. This is used for velocity-step test and similar rotary tests.

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The **Position p** can be set between $-180^\circ \dots 180^\circ$. **Velocity v** sets the velocity of the movement. **Acceleration a** sets the acceleration until **Velocity v** is reached and the deceleration from **Velocity v** until the chair stops.

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The **Velocity v** can be set between $-200^\circ/s \dots 200^\circ/s$. If the velocity is negative, the chair moves anti-clockwise, while at positive velocity-values it moves clockwise. **Acceleration a** sets the acceleration until this velocity is reached and the deceleration from the **Velocity v** until the chair stops.

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