

Interacoustics

Software Suites



Increase patient satisfaction and clinic efficiency



Interacoustics[®]

leading diagnostic solutions


Interacoustics

Software Suites

The Interacoustics suite concept

Interacoustics provides one of the most comprehensive software based solutions for audiologic diagnostics and hearing aid fitting. The Interacoustics Suites integrate various

measurement modules in a common yet intuitive and user friendly way to enable viewing, storing and printing of shared data.



Diagnostic Suite: Supports all standalone, non-PC based audiometers and tympanometers.

Titan Suite: Supports the Titan Middle Ear Analyzer.

Affinity^{2.0}/Equinox^{2.0} Suite: Supports PC based Audiometry, Hearing Instrument Test and Real Ear with Visible Speech Mapping.


Additional modules: Simulators, speech from harddrive, TEN(HL) test, high frequency, etc.

- *Seamless integration*
- *Easy workflow*
- *Intuitive interfaces*
- *Future safe*


Unique Print Wizard for Combined Reports

Interacoustics Software Suites include a common embedded print wizard. This is a flexible tool that supports individual and clinic reporting requirements to streamline clinic workflow. The Print Wizard combines audiometric, impedance and

fitting data for full patient reporting. Unlimited customized templates permit useful alternative templates for different test applications. Printing the reports to PDF allows useful integration into most EMR programs.



Merge and configure input from different sources in the Print Wizard



Insert bitmaps such as clinic logo

Individual session dates

Retrieve patient demographics from OtoAccess™ and NOAH

Merge audiometric and impedance test data on a single page

Text field for session observations and conclusions

Highly configurable symbols for quick reference

EMR and 3rd Party Integration

Interacoustics Suites are prepared for the demands of modern Electronic Medical Records and 3rd party systems by using well-defined and documented interfaces with standard protocols like XML and GDT formats. New methods will continually be implemented as new trends arise.

Future Safe Platform

Interacoustics Suites are future proof platforms serving as the core for all future software based products. Users will benefit from this concept for years to come.

Data Storage

The Interacoustics Suites are designed to integrate data through NOAH or our own proprietary OtoAccess database. In addition, all instruments may be run in 'stand-alone' mode and then integrating data via PDF to EMR.

The Interacoustics software suites

Supported Instruments and Modules

Suite	Modules	Instruments
Diagnostic Suite	Audiometry Tympanometry	AS608/AS608e AD226 AD229b/e AC33 AC40 MT10 AT235/AT235h (old/new)
Titan Suite	Tympanometry	Titan
Affinity ^{2.0} Suite	Audiometry REM HIT	Affinity ^{2.0}
Equinox ^{2.0} Suite	Audiometry	Equinox ^{2.0}

Additional optional modules and simulators for the Affinity^{2.0} and Equinox^{2.0} Suites: SISI, Speech from hard drive, High-frequency audiometry (up to 20kHz), Master hearing aid, Hearing loss simulation built, Multi-frequency and TEN(HL) test.

Diagnostic Suite is Windows® 7 and XP compatible

Interacoustics – the best choice

With over 40 years of experience, Interacoustics is dedicated to supplying its customers with the best possible solutions for their audiologic needs. This is accomplished by maintaining a continuous dialogue with healthcare professionals working in all sectors of audiology. Our equipment meets the highest possible engineering standards and we provide design know-how that can only come from close contact with clinical practice.

Solutions on every scale

Designing equipment for every size of clinic around the world puts us in a unique position to offer solutions that fit your exact requirements. Everything from audiometry, tympanometry, electrophysiology, hearing aid testing, balance investigation to rehab are all within our scope of practice.

Designed for diagnostics

We design diagnostic equipment to make interpretations and ultimately – your life easier. This means superior interfaces, well designed screen layouts, professional reports and interaction through networked databases and EMR systems.

Worldwide Support

The Interacoustics name is not only your guarantee for quality and functionality, but also for support. We operate in over 100 countries worldwide through a well coordinated network of distributors and service centers to ensure that you receive superior support and service.

Sales and service in your area:

Interacoustics A/S

Phone: +45 6371 3555 · Fax: +45 6371 3522
E-mail: info@interacoustics.com
DK-5610 Assens, Denmark
Web: www.interacoustics.com

Interacoustics

- Full line of hearing and balance:

- Audiometry
- Combination audiometry & impedance
- Impedance
- ABR
- ASSR
- OAE
- VEMP
- VNG
- Fitting system
- Database

On our website you can find a full description of all our products including technical specifications and product comparisons.



1 - 05/2010

Read more here:
www.interacoustics.com



Interacoustics®

leading diagnostic solutions